


THRB Antibody (F53852)

Catalog No.	Formulation	Size
F53852-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P10828
Applications	Western Blot : 1:1000-2000
Limitations	This THRB antibody is available for research use only.



Western blot testing of human 1) liver, 2) HeLa, 3) PC-3, 4) SK-BR-3 and 5) TT cell lysate with THRB antibody at 1:2000. Predicted molecular weight: 52 kDa.

Description

Thyroid hormone receptor beta is a nuclear hormone receptor that can act as a repressor or activator of transcription. High affinity receptor for thyroid hormones, including triiodothyronine and thyroxine.

Application Notes

Titration of the THRB antibody may be required due to differences in protocols and substrate sensitivity.

Immunogen

A portion of amino acids 23-55 from the human protein was used as the immunogen for the THRB antibody.

Storage

Aliquot the THRB antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.