


Thrombospondin Antibody (R33775)

Catalog No.	Formulation	Size
R33775-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7057
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This Thrombospondin antibody is available for research use only.



Western blot testing of HEK293 lysate overexpressing human THBS1-MYC with Thrombospondin antibody [1ug/ml] in Lane A and probed with anti-MYC [1/1000] in lane C. Mock-transfected HEK293 probed with Thrombospondin antibody [1ug/ml] in Lane B.

Description

Additional name(s) for this target protein: THBS1, THBS, TSP, TSP1; THBS1

Application Notes

Optimal dilution of the Thrombospondin antibody should be determined by the researcher.

Immunogen

Amino acids EGEARETKACKKD were used as the immunogen for this Thrombospondin antibody.

Storage

Aliquot and store the Thrombospondin antibody at -20oC.