


# **Thyroid Peroxidase Antibody (Biotin Conjugate) (R36270BTN)**

Catalog No.	Formulation	Size
R36270BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P07202
Gene ID	7173
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Thyroid Peroxidase antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Thyroperoxidase; TPO

#### **Application Notes**

Optimal dilution of the Thyroid Peroxidase antibody should be determined by the researcher.

1. This Thyroid Peroxidase antibody is specific for isoforms a/e.

# Immunogen

Amino acids TRHVIQVSNEVVTDD were used as the immunogen for this Thyroid Peroxidase antibody.

### **Storage**

Aliquot and store the Thyroid Peroxidase antibody at -20oC.