


## TIMP1 Antibody (R34917)

Catalog No.	Formulation	Size
R34917-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7076
Localization	Cytoplasmic, secreted
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TIMP1 antibody is available for research use only.



Western blot testing of human kidney lysate with TIMP1 antibody at 0.1ug/ml. Expected molecular weight: 23-28 kDa depending on the level of glycosylation.

## Description

Additional name(s) for this target protein: TIMP metalloproteinase inhibitor 1

## Application Notes

Optimal dilution of the TIMP1 antibody should be determined by the researcher.

## **Immunogen**

Amino acids NSLSLAQRRGFTKTY were used as the immunogen for this TIMP1 antibody.

## **Storage**

Aliquot and store the TIMP1 antibody at -20oC.