


## TIMP1 Antibody (RQ5603)

Catalog No.	Formulation	Size
RQ5603	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit IgG
<b>Purity</b>	Affinity purified
<b>Buffer</b>	Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide
<b>UniProt</b>	P01033
<b>Applications</b>	Western Blot : 0.25-0.5ug/ml Direct ELISA : 0.1-0.5ug/ml
<b>Limitations</b>	This TIMP1 antibody is available for research use only.



Western blot testing of 1) human SK-OV-3, 2) human 22RV1, 3) rat ovary and 4) mouse ovary lysate with TIMP1 antibody. Expected molecular weight: 23-28 kDa depending on the level of glycosylation.

## Description

TIMP metalloproteinase inhibitor 1, also known as TIMP1, a tissue inhibitor of metalloproteinases, is a glycoprotein that is expressed from the several tissues of organisms. This protein is a member of the TIMP family. The glycoprotein is a natural inhibitor of the matrix metalloproteinases (MMPs), a group of peptidases involved in degradation of the extracellular matrix. In addition to its inhibitory role against most of the known MMPs, the encoded protein is able to promote cell proliferation in a wide range of cell types, and may also have an anti-apoptotic function. Transcription of this gene is highly inducible in response to many cytokines and hormones. TIMP was found to be located about 22 cM

proximal to OTC (300461).

## **Application Notes**

Optimal dilution of the TIMP1 antibody should be determined by the researcher.

## **Immunogen**

A human recombinant protein (amino acids H30-S164) was used as the immunogen for the TIMP1 antibody.

## **Storage**

After reconstitution, the TIMP1 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.