


# **TIMP3 Antibody (F55019)**


Catalog No.	Formulation	Size
F55019-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F55019-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

## **Bulk quote request**


Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity purified
UniProt	P35625
Applications	Flow Cytometry: 1:10-1:50 (1x10e6 cells) Immunohistochemistry (FFPE): 1:10-1:50 Western Blot: 1:500-1:1000
Limitations	This TIMP3 antibody is available for research use only.




IHC testing of FFPE human breast carcinoma tissue with TIMP3 antibody. HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.



Western blot testing of human HepG2 cell lysate with TIMP3 antibody. Expected molecular weight: 24-33 kDa depending on glycosylation level.





Flow cytometry testing of human MCF7 cells with TIMP3 antibody; Blue=isotype control, Green= TIMP3 antibody.

### **Description**

TIMP3 is one of inhibitors of the matrix metalloproteinases, a group of peptidases involved in degradation of the extracellular matrix (ECM).

#### **Application Notes**

The stated application concentrations are suggested starting points. Titration of the TIMP3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

#### **Immunogen**

A portion of amino acids 113-139 from the human protein was used as the immunogen for the TIMP3 antibody.

#### **Storage**

Aliquot the TIMP3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.