


TIP30 Antibody (R33127)

Catalog No.	Formulation	Size
R33127-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q9BUP3
Gene ID	10553
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TIP30 antibody is available for research use only.



Western blot testing of HeLa lysate with TIP30 antibody at 0.6ug/ml. Predicted molecular weight: ~27/15/30kDa (isoforms 1/2/3).

Description

Additional name(s) for this target protein: HIV-1 Tat interactive protein 2, 30kDa; HTATIP2

Application Notes

Optimal dilution of the TIP30 antibody should be determined by the researcher.

Immunogen

Amino acids AETEALSKLREDFR were used as the immunogen for this TIP30 antibody. This sequence is common to isoforms 1-3.

Storage

Aliquot and store the TIP30 antibody at -20oC.