


## TIP47 Antibody (R36268)

Catalog No.	Formulation	Size
R36268-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10226
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TIP47 antibody is available for research use only.



Western blot testing of A) HeLa and B) human colon cancer lysate with TIP47 antibody at 0.03ug/ml. Image shows isoforms 1-3 being detected.

## Description

Additional name(s) for this target protein: Tail interacting protein 47 kD; perilipin 3; PLIN3

## Application Notes

Optimal dilution of the TIP47 antibody should be determined by the researcher.

## Immunogen

Amino acids TRGAVQSGVDKTKS were used as the immunogen for this TIP47 antibody. This sequence is common to

isoforms 1-3.

## **Storage**

Aliquot and store the TIP47 antibody at -20oC.