



## Tip60 Antibody / 60 kDa Tat-interactive protein / KAT5 (F40965)

Catalog No.	Formulation	Size
F40965-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40965-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse
<b>Predicted Reactivity</b>	Rat
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Purified
<b>UniProt</b>	Q92993
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This Tip60 antibody is available for research use only.



## Description

HTATIP/Tip60 belongs to the MYST family of histone acetyl transferases (HATs) and was originally isolated as an HIV-1 TAT-interactive protein. HATs play important roles in regulating chromatin remodeling, transcription and other nuclear processes by acetylating histone and nonhistone proteins. This protein is a histone acetylase that has a role in DNA repair and apoptosis and is thought to play an important role in signal transduction.

## Application Notes

Titration of the Tip60 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 33-64 from the human protein was used as the immunogen for this Tip60 antibody.

## Storage

Aliquot the Tip60 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.