


TIPRL Antibody (R35633)

Catalog No.	Formulation	Size
R35633-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	261726
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TIPRL antibody is available for research use only.



Western blot testing of K562 lysate with TIPRL antibody at 1ug/ml. Predicted/observed molecular weight: ~31kDa.

Description

Additional name(s) for this target protein: TOR signaling pathway regulator-like

Application Notes

Optimal dilution of the TIPRL antibody should be determined by the researcher.

1. This TIPRL antibody is specific for isoform 1.

Immunogen

Amino acids RIDPNPADSQKSTQ were used as the immunogen for this TIPRL antibody.

Storage

Aliquot and store the TIPRL antibody at -20oC.