


TIRAP Antibody [isoform b] (R33776)

Catalog No.	Formulation	Size
R33776-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	114609
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TIRAP antibody is available for research use only.



Western blot testing of human liver lysate with TIRAP antibody at 1ug/ml. Predicted molecular weight: ~25kDa, observed here at ~30kDa.

Description

Additional name(s) for this target protein: Toll-interleukin 1 receptor (TIR) domain-containing adapter protein; Mal

Application Notes

Optimal dilution of the TIRAP antibody should be determined by the researcher.

1. This TIRAP antibody is specific to isoform b of the protein.

Immunogen

Amino acids EGEGERDSATVSDL were used as the immunogen for this TIRAP antibody. This sequence is specific for isoform b.

Storage

Aliquot and store the TIRAP antibody at -20oC.