


TIRAP Antibody (R35486)

Catalog No.	Formulation	Size
R35486-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	114609
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TIRAP antibody is available for research use only.



Western blot testing of K562 cell lysate with TIRAP antibody at 0.1ug/ml. Predicted molecular weight: ~24/25kDa (isoforms a/b).

Description

Additional name(s) for this target protein: Toll-interleukin 1 receptor (TIR) domain containing adaptor protein

Application Notes

Optimal dilution of the TIRAP antibody should be determined by the researcher.

Immunogen

Amino acids QTLLKKPKKRPNSPE were used as the immunogen for this TIRAP antibody. This sequence is common to

isoforms a and b.

Storage

Aliquot and store the TIRAP antibody at -20oC.