


TLE4 Antibody [clone 384CT24.1.3.4] (F53481)

Catalog No.	Formulation	Size
F53481-0.1ML	In ascites with 0.09% sodium azide	0.1 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Ascites
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgM
Clone Name	384CT24.1.3.4
Purity	Ascites
UniProt	Q04727
Applications	Western Blot : 1:100-1:3200
Limitations	This TLE4 antibody is available for research use only.



Description

Transcriptional corepressor that binds to a number of transcription factors. Inhibits the transcriptional activation mediated by PAX5, and by CTNNB1 and TCF family members in Wnt signaling. The effects of full-length TLE family members may be modulated by association with dominant-negative AES (By similarity).

Application Notes

Titration of the TLE4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 262-290 from the human protein was used as the immunogen for this TLE4 antibody.

Storage

Aliquot the TLE4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.