


TMPRSS3 Antibody (R33830)

Catalog No.	Formulation	Size
R33830-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64699
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TMPRSS3 antibody is available for research use only.



Western blot testing of A431 cell lysate with TMPRSS3 antibody at 0.1ug/ml. Predicted molecular weight: ~49kDa.

Description

Additional name(s) for this target protein: Transmembrane protease, serine 3

Application Notes

Optimal dilution of the TMPRSS3 antibody should be determined by the researcher.

Immunogen

Amino acids EKIVYHSKYKPKR were used as the immunogen for this TMPRSS3 antibody.

Storage

Aliquot and store the TMPRSS3 antibody at -20oC.