


TMPRSS4 Antibody (R33647)

Catalog No.	Formulation	Size
R33647-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	56649
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This TMPRSS4 antibody is available for research use only.



Western blot testing of human bone marrow lysate with TMPRSS4 antibody at 1ug/ml.
Predicted molecular weight: ~48kDa.

Description

Additional name(s) for this target protein: Transmembrane protease, serine 4

Application Notes

Optimal dilution of the TMPRSS4 antibody should be determined by the researcher.

Immunogen

Amino acids DSTRCNADDAYQGE were used as the immunogen for this TMPRSS4 antibody.

Storage

Aliquot and store the TMPRSS4 antibody at -20oC.