


TMX1 Antibody (R35770)


Catalog No.	Formulation	Size
R35770-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	81542
Applications	Western Blot : 0.01-0.03ug/ml Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This TMX1 antibody is available for research use only.



Western blot testing of human liver lysate with TMX1 antibody at 0.01ug/ml. Predicted molecular weight: ~32kDa.



IHC staining of FFPE human spleen with TMX1 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Thioredoxin domain containing 1, TXNDC1, TXNDC, TMX

Application Notes

Optimal dilution of the TMX1 antibody should be determined by the researcher.

Immunogen

Amino acids RSLGPSLATDKS were used as the immunogen for this TMX1 antibody.

Storage

Aliquot and store the TMX1 antibody at -20oC.