


TNFR1 Antibody (R35851)


Catalog No.	Formulation	Size
R35851-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7132
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TNFR1 antibody is available for research use only.



Western blot testing of human lung lysate with TNFR1 antibody at 0.3ug/ml. Predicted molecular weight ~51kDa.



Western blot testing of A549 (A), HepG2 (B) and K562 (C) lysates with TNFR1 antibody at 0.5ug/ml. Predicted molecular weight ~51kDa.

Description

Additional name(s) for this target protein: Tumor necrosis factor receptor superfamily, member 1A; TNFRSF1A

Application Notes

Optimal dilution of the TNFR1 antibody should be determined by the researcher.

Immunogen

Amino acids SDHEIDRLELQNGR were used as the immunogen for this TNFR1 antibody.

Storage

Aliquot and store the TNFR1 antibody at -20oC.