


TNNI1 Antibody (R35004)

Catalog No.	Formulation	Size
R35004-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7135
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TNNI1 antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with TNNI1 antibody at 0.03ug/ml.
Predicted molecular weight: ~22kDa.

Description

Additional name(s) for this target protein: Troponin I type 1

Application Notes

Optimal dilution of the TNNI1 antibody should be determined by the researcher.

Immunogen

Amino acids PEVERKPKITASRK were used as the immunogen for this TNNI1 antibody.

Storage

Aliquot and store the TNNI1 antibody at -20oC.