


## TNNI2 Antibody (R34039)

Catalog No.	Formulation	Size
R34039-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7136
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This TNNI2 antibody is available for research use only.



Western blot testing of human [A], mouse [B] and rat [C] skeletal muscle lysate with TNNI2 antibody at 1ug/ml. Predicted molecular weight: ~21kDa, observed here at ~25kDa.

## Description

Additional name(s) for this target protein: Troponin I type 2

## Application Notes

Optimal dilution of the TNNI2 antibody should be determined by the researcher.

## Immunogen

Amino acids EVRVQKTSKELED were used as the immunogen for this TNNI2 antibody.

## **Storage**

Aliquot and store the TNNI2 antibody at -20oC.