


## TNNI2 Antibody (R34628)

Catalog No.	Formulation	Size
R34628-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7136
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TNNI2 antibody is available for research use only.



Western blot testing of human [A], mouse [B] and rat [C] skeletal muscle lysate with TNNI2 antibody at 0.03ug/ml. Predicted molecular weight: ~21kDa, observed here at ~25kDa.

## Description

Additional name(s) for this target protein: Troponin I type 2

## Application Notes

Optimal dilution of the TNNI2 antibody should be determined by the researcher.

## Immunogen

Amino acids KRNRAITARRQHLKS were used as the immunogen for this TNNI2 antibody.

## Storage

Aliquot and store the TNNI2 antibody at -20oC.