


## TNNI3 Antibody (R33976)

Catalog No.	Formulation	Size
R33976-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P19429
Gene ID	7137
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TNNI3 antibody is available for research use only.



Western blot testing of human heart lysate with TNNI3 antibody at 1ug/ml. Predicted molecular weight: ~24kDa, can also observed at ~28kDa.

### Description

Additional name(s) for this target protein: Troponin I type 3

### Application Notes

Optimal dilution of the TNNI3 antibody should be determined by the researcher.

## **Immunogen**

Amino acids ERRGEKGRALST were used as the immunogen for this TNNI3 antibody.

## **Storage**

Aliquot and store the TNNI3 antibody at -20oC.