


TP53 Antibody / p53 [clone DDE-20] (RQ5062)

Catalog No.	Formulation	Size
RQ5062	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	DDE-20
Purity	Affinity purified
UniProt	P04637
Applications	Western Blot : 1:500-1:2000
Limitations	This TP53 antibody is available for research use only.



Western blot testing of human 1) Raji and 2) HepG2 cell lysate with TP53 antibody. Predicted molecular weight ~53 kDa.

Description

The TP53 gene encodes a tumor suppressor protein containing transcriptional activation, DNA binding, and oligomerization domains. [RefSeq]

Researchers studying p53 expression and TP53 signaling can explore our Protein Microarray Validated [p53 Antibody](#) page for additional reagents targeting this key tumor suppressor in cancer, DNA damage response and cell cycle research.

Application Notes

Optimal dilution of the TP53 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human p53 / TP53 was used as the immunogen for the TP53 antibody.

Storage

Store the TP53 antibody at -20oC.