


# TP53 Antibody / p53 [clone PCRP-TP53-1F7] (V8500)


| Catalog No.    | Formulation  | Size   |
|----------------|--|--------|
| V8500-100UG    | 0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide | 100 ug |
| V8500-20UG     | 0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide | 20 ug  |
| V8500SAF-100UG | 1 mg/ml in 1X PBS; BSA free, sodium azide free                             | 100 ug |

## **Bulk quote request**


| Availability       | 1-3 business days                                      |
|--------------------|--|
| Species Reactivity | Human  |
| Format             | Purified   |
| Clonality          | Monoclonal (mouse origin)                              |
| Isotype            | Mouse IgG1, kappa                                      |
| Clone Name         | PCRP-TP53-1F7  |
| Purity             | Protein G affinity chromatography                      |
| UniProt            | P04637   |
| Localization       | Nuclear  |
| Applications       | Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT |
| Limitations        | This TP53 antibody is available for research use only. |



Analysis of HuProt(TM) microarray containing more than 19,000 full-length human proteins using TP53 antibody. These results demonstrate the foremost specificity of the PCRP-TP53-1F7 mAb. Z- and S- score: The Z-score represents the strength of a signal that an antibody (in combination with a fluorescently-tagged anti-IgG secondary Ab) produces when binding to a particular protein on the HuProt(TM) array. Z-scores are described in units of standard deviations (SD's) above the mean value of all signals generated on that array. If the targets on the HuProt(TM) are arranged in descending order of the Z-score, the S-score is the difference (also in units of SD's) between the Z-scores. The S-score therefore represents the relative target specificity of an Ab to its intended target.



IHC staining of FFPE human colon carcinoma with TP53 antibody (clone PCRP-TP53-1F7). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.



SDS-PAGE analysis of purified, BSA-free TP53 antibody (clone PCRP-TP53-1F7) as confirmation of integrity and purity.

## **Description**

Recognizes a 53kDa protein, which is identified as p53 suppressor gene product. It reacts with the mutant as well as the wild form of p53 protein. p53 is a tumor suppressor gene expressed in a wide variety of tissue types and is involved in regulating cell growth, replication, and apoptosis. It binds to MDM2, SV40 T antigen and human papilloma virus E6 protein. Positive nuclear staining with p53 antibody has been reported to be a negative prognostic factor in breast carcinoma, lung carcinoma, colorectal, and urothelial carcinoma. Anti-p53 positivity has also been used to differentiate uterine serous carcinoma from endometrioid carcinoma as well as to detect intratubular germ cell neoplasia. Mutations involving p53 are found in a wide variety of malignant tumors, including breast, ovarian, bladder, colon, lung, and melanoma.

## **Application Notes**

Optimal dilution of the TP53 antibody should be determined by the researcher.

#### **Immunogen**

Recombinant human full-length TP53 protein was used as the immunogen for the TP53 antibody.

#### **Storage**

Store the TP53 antibody at 2-8oC (with azide) or aliquot and store at -20oC or colder (without azide).