


TPD52 Antibody (F53908)

Catalog No.	Formulation	Size
F53908-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P55327
Applications	Western Blot : 1:1000-2000
Limitations	This TPD52 antibody is available for research use only.



Western blot testing of 1) A549, 2) HCT116, 3) MCF-7, 4) NCI-H1299, 5) Raji and 6) THP-1 cell lysate with TPD52 antibody at 1:2000. Expected molecular weight: 20~26 kDa (isoforms 1-7).

Description

Alternate names: Tumor protein D52; N8

Application Notes

Titration of the TPD52 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 193-223 from the human protein was used as the immunogen for the TPD52 antibody.

Storage

Aliquot the TPD52 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.