

TPD52L1 Antibody (R33919)

Catalog No.	Formulation	Size
R33919-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7164
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This TPD52L1 antibody is available for research use only.



Western blot testing of human breast cancer lysate with TPD52L1 antibody at 0.1ug/ml.
Predicted molecular weight: ~23/16/15kDa (isoforms 1/3/4).

Description

Additional name(s) for this target protein: Tumor protein D52-like 1

Application Notes

Optimal dilution of the TPD52L1 antibody should be determined by the researcher.

1. This TPD52L1 antibody will detect isoforms 1/3/4 of the protein.

Immunogen

Amino acids EPLQGTDEDAVASAD were used as the immunogen for this TPD52L1 antibody. This sequence is common to isoforms 1, 3 and 4.

Storage

Aliquot and store the TPD52L1 antibody at -20oC.