


## TPH2 Antibody / Tryptophan Hydroxylase 2 (R34844)

Catalog No.	Formulation	Size
R34844-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Dog, Pig, Cow
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	121278
<b>Applications</b>	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:16000
<b>Limitations</b>	This TPH2 antibody is available for research use only.



IHC testing of FFPE human brain stem with TPH2 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

### Description

Additional name(s) for this target protein: Tryptophan hydroxylase 2

### Application Notes

Optimal dilution of the TPH2 antibody should be determined by the researcher.

## **Immunogen**

Amino acids NKPNSGKNDDKGNK were used as the immunogen for this TPH2 antibody.

## **Storage**

Aliquot and store the TPH2 antibody at -20°C.