


TPH2 Antibody / Tryptophan Hydroxylase 2 (R35616)

Catalog No.	Formulation	Size
R35616-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	121278
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TPH2 antibody is available for research use only.



Western blot testing of HEK293 lysate with TPH2 antibody at 0.5ug/ml. Predicted molecular weight: ~56kDa.

Description

Additional name(s) for this target protein: Tryptophan hydroxylase 2

Application Notes

Optimal dilution of the TPH2 antibody should be determined by the researcher.

Immunogen

Amino acids RGFSLDSAVPEEHQ were used as the immunogen for this TPH2 antibody.

Storage

Aliquot and store the TPH2 antibody at -20°C.