


## TRAK1 Antibody (R33234)

Catalog No.	Formulation	Size
R33234-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	22906
Applications	Western Blot : 0.3-1.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This TRAK1 antibody is available for research use only.



Western blot testing of human heart lysate with TRAK1 antibody at 0.3ug/ml. Predicted molecular weight: ~77kDa, observed here at ~90kDa.

### Description

Additional name(s) for this target protein: Trafficking protein, kinesin binding 1

### Application Notes

Optimal dilution of the TRAK1 antibody should be determined by the researcher.

1. This TRAK1 antibody is specific for isoform 1 only.

## **Immunogen**

Amino acids CGAKLSKQTSLR were used as the immunogen for this TRAK1 antibody.

## **Storage**

Aliquot and store the TRAK1 antibody at -20oC.