

## Transglutaminase 2 Antibody (R33953)

Catalog No.	Formulation	Size
R33953-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7052
Applications	Western Blot : 2-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Transglutaminase 2 antibody is available for research use only.



Western blot testing of A) human placenta and B) mouse kidney tissue lysate with Transglutaminase 2 antibody. Predicted molecular weight: ~78 kDa.

### Description

Additional name(s) for this target protein: TGM2; TG2

### Application Notes

Optimal dilution of the Transglutaminase 2 antibody should be determined by the researcher.

### Immunogen

Amino acids ERDLYLENPEIKIR were used as the immunogen for this Transglutaminase 2 antibody. This sequence is specific to isoform a.

## Storage

Aliquot and store the Transglutaminase 2 antibody at -20oC.