


## Transglutaminase 3 Antibody (R36205)

Catalog No.	Formulation	Size
R36205-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7053
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Transglutaminase 3 antibody is available for research use only.



Western blot testing of human tonsil lysate with Transglutaminase 3 antibody at 0.3ug/ml. Predicted/observed molecular weight: ~77kDa.

### Description

Additional name(s) for this target protein: TGM3, TG3

### Application Notes

Optimal dilution of the Transglutaminase 3 antibody should be determined by the researcher.

### Immunogen

Amino acids TLEVLNEARVRKP were used as the immunogen for this Transglutaminase 3 antibody.

## Storage

Aliquot and store the Transglutaminase 3 antibody at -20oC.