


Transglutaminase 6 Antibody (R33519)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R33519-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|-----------------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human, Mouse, Rat |
| Predicted Reactivity | Cow, Dog, Pig |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 343641 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000 |
| Limitations | This Transglutaminase 6 antibody is available for research use only. |



Western blot testing of rat lung lysate with Transglutaminase 6 antibody at 2ug/ml.
Predicted molecular weight: ~80kDa.

Description

Additional name(s) for this target protein: Transglutaminase 6; Tgm6; Tg6

Application Notes

Optimal dilution of the Transglutaminase 6 antibody should be determined by the researcher.

Immunogen

Amino acids DKYVDSYGRTLED were used as the immunogen for this Transglutaminase 6 antibody.

Storage

Aliquot and store the Transglutaminase 6 antibody at -20oC.