


Transmembrane protease serine 2 Antibody / TMPRSS2 [clone TMPRSS2/7420] (V4849)


Catalog No.	Formulation	Size
V4849-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V4849-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V4849SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, kappa
Clone Name	TMPRSS2/7420
Purity	Protein A/G affinity
UniProt	O15393
Localization	Secreted, Cell Surface
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml for 30 min at RT
Limitations	This Transmembrane protease serine 2 antibody is available for research use only.



IHC staining of FFPE human prostate tissue with Transmembrane protease serine 2 antibody (clone TMPRSS2/7420). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.



SDS-PAGE analysis of purified, BSA-free Transmembrane protease serine 2 antibody (clone TMPRSS2/7420) as confirmation of integrity and purity.

Description

Extracellular proteases mediate the digestion of neighboring extracellular matrix components in initial tumor growth, allow desquamation of tumor cells into the surrounding environment, provide the basis for invasion of basement membranes in targeted metastatic organs and are required for release and activation of many growth and angiogenic factors. The TMPRSS2 gene encodes a 492 amino acid multimeric serine protease, which is mainly expressed in the mouse prostate and kidney, and is also expressed in the human small intestine, prostate, colon, stomach and salivary gland. TMPRSS2 contains several domains, including a serine protease domain of the S1 family, a scavenger receptor cysteine-rich domain of group A, an LDL receptor class A domain and a transmembrane domain. TMPRSS2 is expressed as a full length form and a cleaved protease domain and its expression is increased by androgenic hormones. TMPRSS2 is also expressed in prostate carcinoma, suggesting that it may play a role in prostate carcinogenesis.

Application Notes

Optimal dilution of the Transmembrane protease serine 2 antibody should be determined by the researcher.

Immunogen

Recombinant full-length human protein was used as the immunogen for the Transmembrane protease serine 2 antibody.

Storage

Aliquot the Transmembrane protease serine 2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.