


Transthyretin Antibody / TTR / Prealbumin (R36383)


Catalog No.	Formulation	Size
R36383-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7276
Localization	Cytoplasmic, membrane
Applications	IHC (FFPE) : 10ug/ml Enzyme ImmunoAssay : binds recombinant protein (use at 5-10ug/ml) ELISA (peptide) LOD : 1:1000
Limitations	This Transthyretin antibody is available for research use only.



Sandwich ELISA with Transthyretin antibody at 5ug/ml as the detect



IHC testing of FFPE human pancreas tissue with Transthyretin antibody at 10ug/ml. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.

Description

Additional name(s) for this target protein: TTR; prealbumin

Application Notes

Optimal dilution of the Transthyretin antibody should be determined by the researcher.

Immunogen

Amino acids YKVEIDTKSYWK were used as the immunogen for this Transthyretin antibody.

Storage

Aliquot and store the Transthyretin antibody at -20oC.