


## TREML1 Antibody (R33240)

Catalog No.	Formulation	Size
R33240-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	340205
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This TREML1 antibody is available for research use only.



Western blot testing of human bone marrow lysate with TREML1 antibody at 0.5ug/ml.  
Predicted molecular weight: ~33kDa.

## Description

Additional name(s) for this target protein: Triggering receptor expressed on myeloid cells-like 1

## Application Notes

Optimal dilution of the TREML1 antibody should be determined by the researcher.

## Immunogen

Amino acids CGPAQNPPNNQTPS were used as the immunogen for this TREML1 antibody.

## **Storage**

Aliquot and store the TREML1 antibody at -20oC.