


TREX2 Antibody (R34279)

Catalog No.	Formulation	Size
R34279-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11219
Localization	Nuclear
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TREX2 antibody is available for research use only.



IHC testing of FFPE human breast with TREX2 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Three prime repair exonuclease 2

Application Notes

Optimal dilution of the TREX2 antibody should be determined by the researcher.

Immunogen

Amino acids HRSSLENPEHDES were used as the immunogen for this TREX2 antibody.

Storage

Aliquot and store the TREX2 antibody at -20oC.