


## Trib1 Antibody (R34992)

Catalog No.	Formulation	Size
R34992-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10221
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Trib1 antibody is available for research use only.



Western blot testing of NIH3T3 lysate with Trib1 antibody at 0.1ug/ml. Predicted molecular weight: ~41kDa.

## Description

Additional name(s) for this target protein: Tribbles 1

## Application Notes

Optimal dilution of the Trib1 antibody should be determined by the researcher.

## **Immunogen**

Amino acids PEHVSPKARCLIRS were used as the immunogen for this Trib1 antibody.

## **Storage**

Aliquot and store the Trib1 antibody at -20oC.