


TRIM11 Antibody (R33330)

Catalog No.	Formulation	Size
R33330-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	81559
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TRIM11 antibody is available for research use only.



IHC testing of FFPE human kidney with TRIM11 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Tripartite motif-containing 11

Application Notes

Optimal dilution of the TRIM11 antibody should be determined by the researcher.

Immunogen

Amino acids CRPKGGSDDLAPQ were used as the immunogen for this TRIM11 antibody.

Storage

Aliquot and store the TRIM11 antibody at -20oC.