


## TRIM11 Antibody (R33382)

Catalog No.	Formulation	Size
R33382-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	81559
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TRIM11 antibody is available for research use only.



Western blot testing of human spleen lysate with TRIM11 antibody at 0.3ug/ml.  
Predicted molecular weight: ~53kDa.

### Description

Additional name(s) for this target protein: Tripartite motif containing 11

### Application Notes

Optimal dilution of the TRIM11 antibody should be determined by the researcher.

## Immunogen

Amino acids DALRRVQDVKLQPPE were used as the immunogen for this TRIM11 antibody.

## Storage

Aliquot and store the TRIM11 antibody at -20oC.