


TRIM2 Antibody (R36082)

Catalog No.	Formulation	Size
R36082-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23321
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TRIM2 antibody is available for research use only.



Western blot testing of mouse brain lysate with TRIM2 antibody at 0.1ug/ml.
Predicted/observed molecular weight (mouse): ~81kDa.

Description

Additional name(s) for this target protein: Tripartite motif-containing 2

Application Notes

Optimal dilution of the TRIM2 antibody should be determined by the researcher.

Immunogen

Amino acids STFDELQKTLNVRKS were used as the immunogen for this TRIM2 antibody.

Storage

Aliquot and store the TRIM2 antibody at -20oC.