


## TRIM21 Antibody (R33302)

Catalog No.	Formulation	Size
R33302-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6737
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TRIM21 antibody is available for research use only.



Western blot testing of HEK293 lysate overexpressing human TRIM21 with MYC tag with TRIM21 antibody [1ug/ml] in Lane A and probed with anti-MYC tag [1/1000] in lane C. Mock-transfected HEK293 probed with TRIM21 antibody [1ug/ml] in Lane B.

## Description

Additional name(s) for this target protein: Tripartite motif-containing 21

## Application Notes

Optimal dilution of the TRIM21 antibody should be determined by the researcher.

## Immunogen

Amino acids CPLNIGSQGSTDY were used as the immunogen for this TRIM21 antibody.

## **Storage**

Aliquot and store the TRIM21 antibody at -20oC.