


TRIM35 Antibody (R33980)

Catalog No.	Formulation	Size
R33980-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23087
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TRIM35 antibody is available for research use only.



Western blot testing of K562 lysate with TRIM35 antibody at 1ug/ml. Predicted molecular weight: ~56kDa.

Description

Additional name(s) for this target protein: Tripartite motif-containing 35; MAIR

Application Notes

Optimal dilution of the TRIM35 antibody should be determined by the researcher.

Immunogen

Amino acids ERSPDVSPGPSRS were used as the immunogen for this TRIM35 antibody.

Storage

Aliquot and store the TRIM35 antibody at -20oC.