


## TRIM5 Antibody (R33171)

Catalog No.	Formulation	Size
R33171-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	85363
Applications	Western Blot : 2-6ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TRIM5 antibody is available for research use only.



Western blot testing of human placenta lysate with TRIM5 antibody at 2ug/ml. Predicted molecular weight: ~56 kDa (isoform alpha), ~40 kDa (isoform gamma), ~38 kDa (isoform delta), ~31 kDa (isoform epsilon), ~30 kDa (isoform iota).

## Description

Additional name(s) for this target protein: Tripartite motif-containing 5

## Application Notes

Optimal dilution of the TRIM5 antibody should be determined by the researcher.

## Immunogen

Amino acids ASGILVNVKEEVTC were used as the immunogen for this TRIM5 antibody.

## **Storage**

Aliquot and store the TRIM5 antibody at -20oC.