


TRIM5 Antibody (R33290)

Catalog No.	Formulation	Size
R33290-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	85363
Localization	Cytoplasmic, nuclear
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TRIM5 antibody is available for research use only.



IHC testing of FFPE human kidney with TRIM5 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Tripartite motif-containing 5

Application Notes

Optimal dilution of the TRIM5 antibody should be determined by the researcher.

1. This TRIM5 antibody is specific to the alpha isoform.

Immunogen

Amino acids CNIEKNENYQPK were used as the immunogen for this TRIM5 antibody.

Storage

Aliquot and store the TRIM5 antibody at -20oC.