


TRIM71 Antibody (R34841)

Catalog No.	Formulation	Size
R34841-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	131405
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TRIM71 antibody is available for research use only.



Western blot testing of NIH3T3 lysate with TRIM71 antibody at 0.5ug/ml. Predicted molecular weight: ~93kDa. Both observed bands are blocked by addition of the immunizing peptide.

Description

Additional name(s) for this target protein: Tripartite motif-containing 71 antibody, LIN-41, LIN41

Application Notes

Optimal dilution of the TRIM71 antibody should be determined by the researcher.

Immunogen

Amino acids NKLESTISAVQ were used as the immunogen for this TRIM71 antibody.

Storage

Aliquot and store the TRIM71 antibody at -20oC.