


TRIM8 Antibody (R36370)

Catalog No.	Formulation	Size
R36370-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	81603
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TRIM8 antibody is available for research use only.



IHC testing of FFPE human kidney with TRIM8 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Tripartite motif-containing 8; GERP

Application Notes

Optimal dilution of the TRIM8 antibody should be determined by the researcher.

Immunogen

Amino acids YGQPSTKHVTS were used as the immunogen for this TRIM8 antibody.

Storage

Aliquot and store the TRIM8 antibody at -20oC.