


TRIP6 Antibody (R35915)

Catalog No.	Formulation	Size
R35915-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7205
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TRIP6 antibody is available for research use only.



Western blot testing of K562 lysate with TRIP6 antibody at 2ug/ml. Predicted molecular weight: ~50kDa.

Description

Additional name(s) for this target protein: Thyroid hormone receptor interactor 6

Application Notes

Optimal dilution of the TRIP6 antibody should be determined by the researcher.

Immunogen

Amino acids SHGVLQHTQGLP were used as the immunogen for this TRIP6 antibody.

Storage

Aliquot and store the TRIP6 antibody at -20oC.