


TrkA/B/C Antibody / Pan TRK [clone HHH-14] (RQ5057)

Catalog No.	Formulation	Size
RQ5057	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	HHH-14
Purity	Affinity purified
UniProt	P04629, Q16620, Q16288
Applications	Western Blot : 1:1000-1:2000
Limitations	This TrkA/B/C antibody is available for research use only.



Western blot testing of human fetal brain lysate with TrkA/B/C antibody. Expected molecular weight: 85~140 kDa depending on glycosylation level.

Description

TrkA, TrkB and TrkC are members of the neurotrophic tyrosine kinase receptor (NTKR) family.

Application Notes

Optimal dilution of the TrkA/B/C antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human TrkA/B/C was used as the immunogen for the TrkA/B/C antibody.

Storage

Store the TrkA/B/C antibody at -20°C.