


TRPC4 Antibody (R34454)

Catalog No.	Formulation	Size
R34454-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7223
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TRPC4 antibody is available for research use only.



Western blot testing of rat brain lysate with TRPC4 antibody at 0.3ug/ml.
Predicted/observed molecular weight: ~112kDa.

Description

Additional name(s) for this target protein: Transient receptor potential cation channel, subfamily C, member 4

Application Notes

Optimal dilution of the TRPC4 antibody should be determined by the researcher.

Immunogen

Amino acids KHAKEEDSSIDYD were used as the immunogen for this TRPC4 antibody.

Storage

Aliquot and store the TRPC4 antibody at -20oC.