


TRPC4AP Antibody (R34529)

Catalog No.	Formulation	Size
R34529-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26133
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This TRPC4AP antibody is available for research use only.



Western blot testing of U937 lysate with TRPC4AP antibody at 1ug/ml. Predicted molecular weight: ~90kDa.

Description

Additional name(s) for this target protein: Transient receptor potential cation channel, subfamily C, member 4

Application Notes

Optimal dilution of the TRPC4AP antibody should be determined by the researcher.

Immunogen

Amino acids KLTQETTYPNTYIFD were used as the immunogen for this TRPC4AP antibody.

Storage

Aliquot and store the TRPC4AP antibody at -20oC.