


# TRPC6 Antibody (R33837)

Catalog No.	Formulation	Size
R33837-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7225
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TRPC6 antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Transient receptor potential cation channel, subfamily C, member 6

## **Application Notes**

Optimal dilution of the TRPC6 antibody should be determined by the researcher.

#### **Immunogen**

Amino acids EKLSMEPNQEETNR were used as the immunogen for this TRPC6 antibody. **Storage** Aliquot and store the TRPC6 antibody at -20oC.